

MARKES
international



Analisi dell'aria con Thermal Desorber

Davide Facciabene
Product Specialist GC & GC-MS



Pescara 24 Novembre 2015

 **artaabruzzo**
agenzia regionale per la tutela dell'ambiente

Thermal Desorption, cos'è...

Sistema di campionamento ed analisi dell'aria, sia in campo aperto (zone industriali, residenziali, discariche) sia in campo chiuso (luoghi di lavoro, abitati, ecc.)



Thermal Desorption, cos'è...

Il Sistema di campionamento consiste di un tubo metallic o di vetro contenente un sorbente, in cui viene fatta passare l'aria da analizzare in maniera forzata o passiva.

In tal modo, il sorbente estrae la parte organica volatile dell'aria.



Thermal Desorption, cos'è...

Successivamente il tubo viene portato al TD che mediante riscaldamento e sotto il flusso di un gas inerte, libera i principi attivi direttamente nel GC-MS




Campi di applicazione...

For sample matrices that cannot be directly introduced to the analyser (GC) or require pre-concentration

Dilute Environmental Vapour Samples



Material Emissions – construction, automotive




Food, Flavour and Fragrance



Defence and Forensic




Thermal Desorption, cos'è...




Thermal Desorption, cos'è...

Sample passes onto the sorbent



Compounds of interest are adsorbed
on the sorbent surface

Thermal Desorption, cos'è...



Lighter gases such as nitrogen,
argon and carbon dioxide pass
through the sorbent

Thermal Desorption, cos'è...



The sorbent tube is now heated in a reversed flow of clean carrier gas (back flushed)

Thermal Desorption, cos'è...

Compounds are released from the sorbent into the flow of carrier gas



It combines preconcentration, desorption/extraction and GC injection into one sensitive and fully automated operation

It is a simple extension of the technique of Gas Chromatography and is a **sample introduction technology** for difficult or real-world samples. Analytes undergo **pre-concentration** from litre samples to μl vapour band on the GC

Thermal Desorption, cos'è...


PROBLEM: Compounds are released SLOWLY from the sorbent tube



Would lead to **very wide** chromatographic peaks and low sensitivity

Thermal Desorption, cos'è...


ANSWER: All modern thermal desorption instruments incorporate a two-stage desorption procedure.



Thermal Desorption, cos'è...

STAGE 2

Rapid transfer of
compounds from cold
trap to GC



- Cold trap heated rapidly (100°C/sec) for sharp chromatographic peaks
- Transfer of analytes through narrow bore transfer line of cold trap for greater volatility range

Quali composti...


Any volatile or semi-volatile organic compounds which meet the following criteria:

- ✓ $\leq n\text{-C}_{40}$, bpt $\leq 525^\circ\text{C}$
- ✓ Compatible with 'standard' GC analysis
- ✓ The sorbent or matrix containing the compounds is compatible with the high temperatures required
- ✗ Inorganic compounds
- ✗ Most permanent gases (CO_2 , Ar, N_2 etc.) and other compounds with very low boiling points (methane, formaldehyde)
- ✗ Compounds bigger than $n\text{-C}_{40}$
- ✗ Compounds which are not compatible with gas chromatography (or which require on-column injection or derivitisation)


EPA TO-17

1	Propylene	22	Cis-1,2-Dichloroethylene	43	Methyl <i>n</i> -butyl ketone
2	Dichlorodifluoromethane	23	Methyl ethyl ketone	44	Dibromochloromethane
3	1,2-Dichlorotetrafluoroethane	24	Ethyl acetate	45	1,2-Dibromoethane
4	Methyl chloride	25	Tetrahydrofuran	46	Chlorobenzene
5	1,2-Dichloroethane	26	Chloroform	47	Xylene
6	1,3-Butadiene	27	1,1,1-Trichloroethane	48	Xylene
7	Vinyl chloride	28	Cyclohexane	49	Xylene
8	Methyl bromide (bromomethane)	29	Carbon tetrachloride	50	Styrene
9	Chloroethane	30	Benzene	51	Tribromomethane
10	Trichlorotrifluoroethane (Freon® 113)	31	<i>n</i> -Heptane	52	1,1,2,2-Tetrachloroethane
11	Ethanol	32	Trichloroethylene	53	1,2,4-Trimethylbenzene
12	1,2,-Dichloroethylene	33	1,2-Dichloropropane	54	1,3,5-Trimethylbenzene
13	1,1,2-Trichlorotrifluoroethane	34	1,4-Dioxane	55	1-Ethyl-4-methyl benzene
14	Acetone	35	Bromodichloromethane	56	Ethylbenzene
15	Carbon disulfide	36	Trans-1,3-dichloropropene	57	1,2-Dichlorobenzene
16	Isopropyl alcohol	37	Methyl isobutyl ketone	58	1,3-Dichlorobenzene
17	Methylene chloride	38	Toluene	59	Chloromethylbenzene (alpha)
18	Tert-butyl methyl ether	39	Cis-1,3-Dichloropropene	60	1,4-Dichlorobenzene
19	<i>n</i> -Hexane	40	Trans-1,2-Dichloroethylene	61	1,2,4-Trichlorobenzene
20	1,1-Dichloroethane	41	1,1,2-Trichloroethane	62	Hexachloro-1,3-butadiene
21	Vinyl acetate	42	Tetrachloroethylene		

Splitless desorption of
'Air toxics' tube loaded
with 1 L of 1 ppb std
GC/MS



IPA in aria ambientale



Compounds	
Naphthalene	1
1-methyl naphthalene	2
2-methyl naphthalene	3
Acenaphthylene	4
Acenaphthene	5
Fluorene	6
Phenanthrene	7
Anthracene	8
Fluoranthene	9
Pyrene	10
Benz[a]anthracene	11
Chrysene	12
Benzo[b]fluoranthene	13
Benzo[k]fluoranthene	14
Benzo[a]pyrene	15
Indeno[1,2,3-cd]pyrene	16
Dibenzo[a,h]anthracene	17
Benzo[ghi]perylene	18

Thermal Desorption, cos'è...

Sorbent name	Volatility range
Tenax TA	C ₇ – C ₃₀
Carbograph 2TD	C ₈ – C ₂₀
Carbograph 1TD	C _{5/6} – C ₁₄
Carbograph 5TD	C _{3/4} – C _{6/7}
SulfiCarb	C ₃ – C ₈
Carboxen 1003	C ₂ – C ₅
Carbosieve SIII	C ₂ – C ₅

Weak Retention

Water retention

Strong Retention

Campionamento...


Three sampling methods

- a. Passive (diffusive) sampling
- b. Direct sampling
- c. Active (pumped) sampling




Campionamento Passivo (diffusivo)

- Diffusive sampling = a simple and cost effective method of collecting the large number of samples required in many air monitoring programmes. See Markes Application note #8 for an overview
- Vapours migrate across the air gap at a constant “uptake rate” as tube dimensions are consistent (Fick’s law)
- Diffusive sampling is a slow process, typically sample for days



Ambient conc. of vapours in the environment

Zero conc. of vapours at the sorbent surface



Campionamento Passivo (diffusivo)

Only compatible with single bed sorbent tubes - only one end of the tube is exposed

$$\text{Conc (ppm)} = \frac{\text{Mass on tube (ng)}}{\text{Uptake Rate x Sample Time (mins)}}$$

- Uptake rates of many analytes on a range of sorbents have now been published (see Markes Application Note #1 and #42)
- If an uptake rate is not available in the literature it is possible to determine it experimentally but this is difficult
 - Involves collecting several diffusive + pumped samples from the same 'atmosphere'
 - The pumped samples provide you with the known concentration so you can then calculate the uptake rate for the diffusive samplers

Campionamento Passivo (diffusivo)


- Two types of diffusive samplers:
 - Tube-type axial samplers (sorbent tubes)



- Uptake rates typically $1-3\text{ng ppm}^{-1} \text{min}^{-1}$


- Radial Diffusive samplers (sorbent cartridges)

- Use for analysis of low concentrations over shorter periods e.g. workplace monitoring (<1ppb over 8hours)
- Cylindrical sorbent cartridge placed inside empty sorbent tube and analysed
- Uptake rates typically 50-100x faster than axial sampling



Campionamento Attivo (forzato)

- Pump air through sorbent tube
- Flow Rate = 20 – 200 ml/min
- Much faster technique compared to diffusive sampling
- Fully quantitative – know how much is sampled
- **Important** do not exceed breakthrough volume for a compound on a given sorbent



Campionamento Attivo

Constant flow pump

- Pump varies its speed to maintain a constant (programmed) flow
- Set a defined sampling flow and time



Constant pressure pump

- Flow will vary with differing tube impedance
- Must measure / verify flow when changing tube types as different sorbents have differing impedances




Syringe pump

- Simple way to take pumped samples that doesn't require electrical power or re-calibration




Breakthrough



Affected by:

- Type of sorbent **stronger sorbent = stronger interaction**
- Sample volume **lower volume = less risk of breakthrough**
- Temperature **lower temp. = stronger interaction**
- Mass or sorbent **more sorbent = more surface area**



MARKES
international



Grazie per l'attenzione

